

## ProCare Tex 10 W - 20 l

Detergent for wool, liquid concentrate, mildly acidic, 20 l for washing woollens in washing machines.

EAN: 4002516732440 / Material number: 12338960 / Old Material Number: 59706301V1



- No felting of wool fabrics
- Prevents the transfer of dye
- Prevents greying caused by redeposition of soil
- Dermatologically tested

<b>Product classification</b>	
Washing machines	•
<b>Dispensing</b>	
pH value of solution	4,8-5,8
Suitable for automatic dispensing	•
Can be used with all degrees of water hardness	•
<b>Use</b>	
Launderettes	•
Hotels, catering	•
<b>Product properties</b>	
Biodegradable	•
Particularly effective at low temperatures	•
High-level detergency combined with sparing use	•
Density (at 20 °C) in g/cm <sup>3</sup>	1.03
Free from micro-plastics	•
Simple handling	•
<b>Storage</b>	
Minimum storage temperature in °C	0
Maximum storage temperature in °C	40
Store in a dry location that is free of frost	•
Containers should always be re-sealed immediately once the required amount of product has been removed	•
Minimum shelf life in months	24
<b>Dimensions and weight</b>	
External dimensions, net height in mm	385
External dimensions, net width in mm	258
External dimensions, net depth in mm	294
Fill volume in ml	20000
<b>Compliance, certs/symbols</b>	
GHS/CLP compliant	•