

KWT 6722 iS-1

Built-in wine conditioning unit
with FlexiFrame and SommelierSet for demanding wine connoisseurs.

EAN: 4002516712367 / Material number: 12364320 / Old Material Number: 36672202EU1




- Flexible storage of larger wine bottles thanks to FlexiFrame
- Simultaneous wine storage in separate temperature zones
- Odourless refrigeration section thanks to Active AirClean filters
- Present your wine in the perfect light – LED lighting




KWT 6722 iS-1

Built-in wine conditioning unit

with FlexiFrame and SommelierSet for demanding wine connoisseurs.



KWT6722iS, KWT6722iS-1, Installation draw



KWT6722iS, KWT6722iS-1, Interior dimensions (Installation drawings)