

APLW 880
 Retrofitting kit
 for Rec in, fourth water inlet, including booster pump for PLW 8636.



- Required in order to pump a medium out of a buffer container
- For use with an existing on-site recycling tank
- Can be used with existing on-site container, e.g. IBC tank
- Installed by Customer Service

EAN: 4002516745532 / Material number: 12426490 /
 Old Material Number: 69688000D

Product category	
Washer-disinfectors, laboratory	•
Cleaning machines for industrial parts cleaning	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
Product category	
Conversion kit	•
Dimensions and weight	
External dimensions, gross height in mm	123
External dimensions, gross width in mm	315
External dimensions, gross depth in mm	505
Net weight in kg	4.36
Gross weight in kg	4.84