

APLW 883

Retrofitting kit

extra disp. membrane pump & suction lance, int. 5–10 l canisters for PLW 8636.



- For 5–10 l canisters in the dispensing drawer
- The suction lance features level detection
- Installed by Customer Service

EAN: 4002516744603 / Material number: 12426680 /  
Old Material Number: 69688300D

Product category	
Washer-disinfectors, laboratory	•
Cleaning machines for industrial parts cleaning	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
Product category	
Conversion kit	•
Dispenser system	•
Dimensions and weight	
External dimensions, gross height in mm	125
External dimensions, gross width in mm	165
External dimensions, gross depth in mm	350
Net weight in kg	0.79
Gross weight in kg	0.95