

APLW 885

Retrofitting kit

extra disp. membrane pump & suction lance, ext. 10–20 l canisters for PLW 8636.



EAN: 4002516745525 / Material number: 12426770 / Old Material Number: 69688500D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Cleaning machines for industrial parts cleaning	•
<b>Machine affiliation</b>	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
<b>Product category</b>	
Conversion kit	•
Dispenser system	•

- For 10–20 l canisters, positioned outside the machine
- The suction lance features level detection
- Installed by Customer Service