

APLW 886

Retrofitting kit

for dispensing monitoring by oval-wheel counter, for installation in PLW 8636.

EAN: 4002516745433 / Material number: 12426910 / Old Material Number: 69688600D



- More precise quantity measurement using oval-wheel counters
- Installed by Customer Service

Product category	
Washer-disinfectors, laboratory	•
Cleaning machines for industrial parts cleaning	•
Machine affiliation	
PLW 8636 [IPC]	•
PLW 8636 [LAB / LAB MON]	•
Product category	
Conversion kit	•
Dispenser system	•
Dimensions and weight	
External dimensions, gross height in mm	80
External dimensions, gross width in mm	95
External dimensions, gross depth in mm	110
Net weight in kg	0.03
Gross weight in kg	0.08