

ProCare Tex 10 A - 5 l

Sole detergent, liquid concentrate, alkaline, 5 l phosphate-free, with EU Ecolabel.

EAN: 4002516757931 / Material number: 12441960 / Old Material Number: 59701102V1



- For optimal cleaning of various textiles
- Specially developed for Miele Professional
- Based on sustainable raw materials
- Short activation time and high-level detergency with sparing use
- Excellent protein and grease removal
- Awarded the EU Ecolabel

Product category	
Washing machines	•
Wash performance	
Dirt dissolving power	Very good
Protein dissolving power	Good
Grease dissolving power	Good
Dispensing	
Adjustment of dispensing concentration possible	•
pH value of solution	12.0-13.0
Suitable for automatic dispensing	•
Can be used with all degrees of water hardness	•
How to use	
Launderettes	•
Hotels, catering	•
Cleaning sportswear	•
Cleaning outerwear	•
Product properties	
Medium type	Sole detergent
Alkalinity	Alkaline
Aggregate form	Liquid concentrate
Composition	Alkalines, Surfactants
Biodegradable	•
Phosphate-free	•
Effective at temperatures from 30–90 °C	•
Highly effective, even at 40 °C	•
Short activation time	•
High-level detergency combined with sparing use	•
Density (at 20 °C) in g/cm ³	1.17
Based on 100 % sustainable raw materials	•
Storage information	
Minimum storage temperature in °C	0
Maximum storage temperature in °C	40
Store in a dry location that is free of frost	•
Containers should always be re-sealed immediately once the required amount of product has been removed	•
Minimum shelf life in months	24
Dimensions and weight	
External dimensions, net height in mm	256
External dimensions, net width in mm	148
External dimensions, net depth in mm	193
Weight in g	6250
Fill volume in ml	5000
Compliance, test certificates, and symbols	
GHS/CLP compliant	•
EU eco label	•
Verification numbers	
EU Eco label: DE/039/019	•