

ProCare Tex 10 A - 10 l

Sole detergent, liquid concentrate, alkaline, 10 l
for washing white and colourfast coloured items.



- For use in Miele Professional washing machines
- Based on sustainable raw materials
- Short activation time and high-level detergency with sparing use
- Excellent protein and grease removal
- Awarded the EU Ecolabel

EAN: 4002516757771 / Material number: 12441970 / Old Material Number: 59701202V1

Product category	
Washing machines	•
Wash performance	
Dirt dissolving power	Very good
Protein dissolving power	Good
Grease dissolving power	Good
Dispensing	
Adjustment of dispensing concentration possible	•
Suitable for automatic dispensing	•
Can be used with all degrees of water hardness	•
Dispensing recommendation	
pH value at 20 °C	12,6
How to use	
Launderettes	•
Hotels, catering	•
Cleaning sportswear	•
Cleaning outerwear	•
Cleaning workwear	•
Product properties	
Minimum storage temperature in °C	0
Maximum storage temperature in °C	40
Dimensions and weight	
External dimensions, gross depth in mm	228.00
Product properties	
Medium type	Sole detergent
Alkalinity	Alkaline
Aggregate form	Liquid concentrate
Composition	Alkalines, Surfactants
Biodegradable	•
Phosphate-free	•
Effective at temperatures from 30–90 °C	•
Highly effective, even at 40 °C	•
Short activation time	•
High-level detergency combined with sparing use	•
Density (at 20 °C) in g/cm ³	1.17
Based on 100 % sustainable raw materials	•
Storage information	
Minimum storage temperature in °C	0
Maximum storage temperature in °C	40
Store in a dry location that is free of frost	•
Containers should always be re-sealed immediately once the required amount of product has been removed	•
Minimum shelf life in months	24
Dimensions and weight	
External dimensions, net height in mm	314
External dimensions, net width in mm	202
External dimensions, net depth in mm	230
External dimensions, gross height in mm	306
External dimensions, gross width in mm	204
External dimensions, gross depth in mm	228

ProCare Tex 10 A - 10 l

Sole detergent, liquid concentrate, alkaline, 10 l
for washing white and colourfast coloured items.



EAN: 4002516757771 / Material number: 12441970 / Old Material Number: 59701202V1

Weight in g	12500
Gross weight in kg	12.26 kg
Fill volume in ml	10000
Compliance, test certificates, and symbols	
GHS/CLP compliant	•
EU eco label	•
Verification numbers	
EU Eco label: DE/039/019	•