

## ProCare Tex 10 A - 20 l

Sole detergent, liquid concentrate, alkaline, 20 l phosphate-free, with EU Ecolabel.

EAN: 4002516757795 / Material number: 12441980 / Old Material Number: 59701302V1



- For optimal cleaning of various textiles
- Specially developed for Miele Professional
- Based on sustainable raw materials
- Short activation time and high-level detergency with sparing use
- Excellent protein and grease removal
- Awarded the EU Ecolabel

<b>Product category</b>	
Washing machines	•
<b>Wash performance</b>	
Dirt dissolving power	Very good
Protein dissolving power	Good
Grease dissolving power	Good
<b>Dispensing</b>	
Adjustment of dispensing concentration possible	•
pH value of solution	12.0-13.0
Suitable for automatic dispensing	•
Can be used with all degrees of water hardness	•
<b>How to use</b>	
Launderettes	•
Hotels, catering	•
Cleaning sportswear	•
Cleaning outerwear	•
<b>Product properties</b>	
Medium type	Sole detergent
Alkalinity	Alkaline
Aggregate form	Liquid concentrate
Composition	Alkalines, Surfactants
Biodegradable	•
Phosphate-free	•
Effective at temperatures from 30–90 °C	•
Highly effective, even at 40 °C	•
Short activation time	•
High-level detergency combined with sparing use	•
Density (at 20 °C) in g/cm <sup>3</sup>	1.17
Based on 100 % sustainable raw materials	•
<b>Storage information</b>	
Minimum storage temperature in °C	0
Maximum storage temperature in °C	40
Store in a dry location that is free of frost	•
Containers should always be re-sealed immediately once the required amount of product has been removed	•
Minimum shelf life in months	24
<b>Dimensions and weight</b>	
External dimensions, net height in mm	385
External dimensions, net width in mm	258
External dimensions, net depth in mm	294
Weight in g	24000
Fill volume in ml	20000
<b>Compliance, test certificates, and symbols</b>	
GHS/CLP compliant	•
EU eco label	•
<b>Verification numbers</b>	
EU Eco label: DE/039/019	•