

ProCare Tex 10 MA - 5 l

Detergent for delicates/colours, liquid concentrate, mildly alkaline, 5 l highly effective, phosphate-free, with EU Ecolabel.



- For optimal cleaning of delicate textiles
- Specially developed for Miele Professional
- Free from optical brighteners and bleaching agents
- Short activation time, foam sensing
- Good protein and grease removal
- Awarded the EU Ecolabel

EAN: 4002516757665 / Material number: 12442010 / Old Material Number: 59702102V1

Product category	
Washing machines	•
Wash performance	
Dirt dissolving power	Good
Protein dissolving power	Good
Grease dissolving power	Good
Dispensing	
Adjustment of dispensing concentration possible	•
pH value of solution	8.5-9.6
Suitable for automatic dispensing	•
Can be used with all degrees of water hardness	•
How to use	
Hotels, catering	•
Cleaning sportswear	•
Cleaning outerwear	•
Cleaning workwear	•
Product properties	
Medium type	Delicates/coloureds detergent
Alkalinity	Mildly alkaline
Aggregate form	Liquid concentrate
Composition	Alkalines, Preservatives, Surfactants
Phosphate-free	•
Effective at temperatures from 30–90 °C	•
Foam sensing	•
Density (at 20 °C) in g/cm ³	1.03
Storage information	
Minimum storage temperature in °C	4
Maximum storage temperature in °C	40
Store in a dry location that is free of frost	•
Containers should always be re-sealed immediately once the required amount of product has been removed	•
Can be used after defrosting without loss of quality	•
Minimum shelf life in months	24
Dimensions and weight	
External dimensions, net height in mm	256
External dimensions, net width in mm	148
External dimensions, net depth in mm	193
Weight in g	5450
Fill volume in ml	5000
Compliance, test certificates, and symbols	
GHS/CLP compliant	•
EU eco label	•
Verification numbers	
EU Eco label: DE/039/020	•