

## ProCare Tex 10 MA - 20 l

Detergent for delicates/colours, liquid concentrate, mildly alkaline, 20 l highly effective, phosphate-free, with EU Ecolabel.

EAN: 4002516757559 / Material number: 12442070 / Old Material Number: 59702302V1



- For optimal cleaning of delicate textiles
- Specially developed for Miele Professional
- Free from optical brighteners and bleaching agents
- Short activation time, foam sensing
- Good protein and grease removal
- Awarded the EU Ecolabel

<b>Product category</b>	
Washing machines	•
<b>Wash performance</b>	
Dirt dissolving power	Good
Protein dissolving power	Good
Grease dissolving power	Good
<b>Dispensing</b>	
Adjustment of dispensing concentration possible	•
pH value of solution	8.5-9.6
Suitable for automatic dispensing	•
Can be used with all degrees of water hardness	•
<b>How to use</b>	
Hotels, catering	•
Cleaning sportswear	•
Cleaning outerwear	•
Cleaning workwear	•
<b>Product properties</b>	
Medium type	Delicates/coloureds detergent
Alkalinity	Mildly alkaline
Aggregate form	Liquid concentrate
Composition	Alkalines, Preservatives, Surfactants
Phosphate-free	•
Effective at temperatures from 30–90 °C	•
Foam sensing	•
Density (at 20 °C) in g/cm <sup>3</sup>	1.03
<b>Storage information</b>	
Minimum storage temperature in °C	4
Maximum storage temperature in °C	40
Store in a dry location that is free of frost	•
Containers should always be re-sealed immediately once the required amount of product has been removed	•
Can be used after defrosting without loss of quality	•
Minimum shelf life in months	24
<b>Dimensions and weight</b>	
External dimensions, net height in mm	385
External dimensions, net width in mm	258
External dimensions, net depth in mm	294
Weight in g	21000
Fill volume in ml	20000
<b>Compliance, test certificates, and symbols</b>	
GHS/CLP compliant	•
EU eco label	•
<b>Verification numbers</b>	
EU Eco label: DE/039/020	•