

ProCare Tex 10 B+ - 5 l
 Booster, liquid concentrate, neutral, 5 l
 for effective removal of grease stains.



- Gentle degreaser for any type of laundry
- Particularly gentle on fabrics and colours
- Full washing power even at 30 °C
- Can also be used as an emulsifier and wetting agent

EAN: 4002516776185 / Material number: 12495130 / Old Material Number: 59709101V1

Product category	
Washing machines	•
Wash performance	
Dirt dissolving power	Very good
Protein dissolving power	Very good
Grease dissolving power	Very good
Degree of whiteness	Very good
Dispensing	
Suitable for automatic dispensing	•
Manual dispensing	•
Can be used with all degrees of water hardness	•
Dispensing recommendation	
pH value at 20 °C	7,5
How to use	
Launderettes	•
Hotels, catering	•
Cleaning sportswear	•
Cleaning outerwear	•
Cleaning workwear	•
Product properties	
Medium type	Booster
Alkalinity	pH-neutral
Aggregate form	Liquid concentrate
Composition	Surfactants
Biodegradable	•
Phosphate-free	•
Effective at temperatures from 30–90 °C	•
Particularly effective at low temperatures	•
Highly effective, even at 40 °C	•
Short activation time	•
Foam sensing	•
High-level detergency combined with sparing use	•
Density (at 20 °C) in g/cm ³	1.00
Free from micro-plastics	•
Simple handling	•
Ideal additive for all detergents	•
Storage information	
Minimum storage temperature in °C	0
Maximum storage temperature in °C	40
Store in a dry location that is free of frost	•
Containers should always be re-sealed immediately once the required amount of product has been removed	•
Can be used after defrosting without loss of quality	•
Minimum shelf life in months	36
Compliance, test certificates, and symbols	
GHS/CLP compliant	•