

## ProCare Tex 11 COL - 8 kg

Detergent for coloured items, powder, mildly alkaline, 8 kg for colour-fast cleaning of coloured items.



- For use in Miele Professional washing machines
- Highly concentrated detergent for coloured items
- Powerful grease and dirt removal that is gentle on materials and colours
- Full washing power even at 30 °C
- Very good protection against dye transfer
- Dermatologically tested

EAN: 4002516775980 / Material number: 12495170 / Old Material Number: 59710201V1

<b>Product category</b>	
Washing machines	•
<b>Wash performance</b>	
Dirt dissolving power	Very good
Protein dissolving power	Very good
Grease dissolving power	Very good
Antistatic action	Very good
<b>Dispensing</b>	
Manual dispensing	•
Can be used with all degrees of water hardness	•
<b>Dispensing recommendation</b>	
pH value at 20 °C	10,5
<b>How to use</b>	
Launderettes	•
Hotels, catering	•
Cleaning sportswear	•
Cleaning outerwear	•
Cleaning workwear	•
<b>Product properties</b>	
Minimum storage temperature in °C	-20
Maximum storage temperature in °C	40
<b>Product properties</b>	
Medium type	Detergent for delicates/coloured items
Alkalinity	Mildly alkaline
Aggregate form	Powder
Composition	Surfactants, Fragrances, Enzymes, Sequestering agent
Biodegradable	•
Biologically degradable thanks to special-purpose enzymes	•
Phosphate-free	•
Particularly effective at low temperatures	•
Highly effective, even at 40 °C	•
Short activation time	•
High-level detergency combined with sparing use	•
Density (at 20 °C) in g/cm <sup>3</sup>	0.8
Degrades without any residues	•
Free from micro-plastics	•
Simple handling	•
Ideal additive for all detergents	•
<b>Storage information</b>	
Minimum storage temperature in °C	-20
Maximum storage temperature in °C	40
Store in a dry location that is free of frost	•
Containers should always be re-sealed immediately once the required amount of product has been removed	•
Can be used after defrosting without loss of quality	•
Minimum shelf life in months	24
<b>Dimensions and weight</b>	
External dimensions, net height in mm	245
External dimensions, net width in mm	290

## ProCare Tex 11 COL - 8 kg

Detergent for coloured items, powder, mildly alkaline, 8 kg  
for colour-fast cleaning of coloured items.

EAN: 4002516775980 / Material number: 12495170 / Old Material Number: 59710201V1



External dimensions, net depth in mm	290
Weight in g	8000
Gross weight in kg	8.40 kg
<b>Compliance, test certificates, and symbols</b>	
GHS/CLP compliant	•