

ProCare Tex 11 COL - 20 kg
Detergent for colours, powder, mildly alkaline, 20 kg
for colour-fast cleaning of coloured items.



- Highly concentrated cleaning agent for coloured items
- Powerful grease and dirt removal that is gentle on materials and colours
- Full washing power even at 30 °C
- Very good protection against dye transfer
- Dermatologically tested

EAN: 4002516776130 / Material number: 12495180 / Old Material Number: 59710301V1

Product category	
Washing machines	•
Wash performance	
Dirt dissolving power	Very good
Protein dissolving power	Very good
Grease dissolving power	Very good
Antistatic action	Very good
Dispensing	
Manual dispensing	•
Can be used with all degrees of water hardness	•
Dispensing recommendation	
pH value at 20 °C	10,5
How to use	
Dry cleaners	•
Launderettes	•
Hotels, catering	•
Cleaning sportswear	•
Cleaning outerwear	•
Cleaning workwear	•
Product properties	
Medium type	Detergent for delicates/coloured items
Alkalinity	Mildly alkaline
Aggregate form	Powder
Composition	Surfactants, Fragrances, Enzymes, Sequestering agent
Biodegradable	•
Biologically degradable thanks to special-purpose enzymes	•
Phosphate-free	•
Particularly effective at low temperatures	•
Highly effective, even at 40 °C	•
Short activation time	•
High-level detergency combined with sparing use	•
Density (at 20 °C) in g/cm ³	0.8
Degrades without any residues	•
Free from micro-plastics	•
Simple handling	•
Ideal additive for all detergents	•
Storage information	
Minimum storage temperature in °C	-20
Maximum storage temperature in °C	40
Store in a dry location that is free of frost	•
Containers should always be re-sealed immediately once the required amount of product has been removed	•
Can be used after defrosting without loss of quality	•
Minimum shelf life in months	24
Compliance, test certificates, and symbols	
GHS/CLP compliant	•