

ProCare Tex 11 COL - 20 kg  
Detergent for colours, powder, mildly alkaline, 20 kg  
for colour-fast cleaning of coloured items.



- For use in Miele Professional washing machines
- Highly concentrated cleaning agent for coloured items
- Powerful grease and dirt removal that is gentle on materials and colours
- Full washing power even at 30 °C
- Very good protection against dye transfer
- Dermatologically tested

EAN: 4002516776130 / Material number: 12495180 / Old Material Number: 59710301V1

<b>Product category</b>	
Washing machines	•
<b>Wash performance</b>	
Dirt dissolving power	Very good
Protein dissolving power	Very good
Grease dissolving power	Very good
Antistatic action	Very good
<b>Dispensing</b>	
Manual dispensing	•
Can be used with all degrees of water hardness	•
<b>Dispensing recommendation</b>	
pH value at 20 °C	10,5
<b>How to use</b>	
Launderettes	•
Hotels, catering	•
Cleaning sportswear	•
Cleaning outerwear	•
Cleaning workwear	•
<b>Product properties</b>	
Minimum storage temperature in °C	-20
Maximum storage temperature in °C	40
<b>Dimensions and weight</b>	
External dimensions, gross depth in mm	620.00
<b>Product properties</b>	
Medium type	Detergent for delicates/ coloured items
Alkalinity	Mildly alkaline
Aggregate form	Powder
Composition	Surfactants, Fragrances, Enz ymes, Sequestering agent
Biodegradable	•
Biologically degradable thanks to special-purpose enzymes	•
Phosphate-free	•
Particularly effective at low temperatures	•
Highly effective, even at 40 °C	•
Short activation time	•
High-level detergency combined with sparing use	•
Density (at 20 °C) in g/cm <sup>3</sup>	0.8
Degrades without any residues	•
Free from micro-plastics	•
Simple handling	•
Ideal additive for all detergents	•
<b>Storage information</b>	
Minimum storage temperature in °C	-20
Maximum storage temperature in °C	40
Store in a dry location that is free of frost	•
Containers should always be re-sealed immediately once the required amount of product has been removed	•
Can be used after defrosting without loss of quality	•
Minimum shelf life in months	24

ProCare Tex 11 COL - 20 kg  
Detergent for colours, powder, mildly alkaline, 20 kg  
for colour-fast cleaning of coloured items.



EAN: 4002516776130 / Material number: 12495180 / Old Material Number: 59710301V1

#### Dimensions and weight

External dimensions, net height in mm	480
External dimensions, net width in mm	330
External dimensions, net depth in mm	130
External dimensions, gross height in mm	140
External dimensions, gross width in mm	440
External dimensions, gross depth in mm	620
Weight in g	20000
Gross weight in kg	20.25 kg

#### Compliance, test certificates, and symbols

GHS/CLP compliant

•