

ProCare Tex 30 C - 5 l

Neutralising agent, liquid concentrate, acidic, 5 l
for optimum protection of textiles thanks to reliable neutralisation.



- For use in Miele Professional washing machines
- Neutralising agent based on citric acid
- For acidic pre-wash or for use after main wash
- Suitable for all wash processes
- Prevents yellow stains
- Stable, neutral pH values in the textiles

EAN: 4002516776215 / Material number: 12495200 / Old Material Number: 59711101V1

Product category	
Washing machines	•
Dispensing	
Suitable for automatic dispensing	•
Manual dispensing	•
Can be used with all degrees of water hardness	•
Dispensing recommendation	
pH value at 20 °C	2,5
How to use	
Hotels, catering	•
Cleaning sportswear	•
Cleaning outerwear	•
Cleaning workwear	•
Product properties	
Minimum storage temperature in °C	0
Maximum storage temperature in °C	40
Dimensions and weight	
External dimensions, gross depth in mm	193.00
Product properties	
Medium type	Neutralising agents
Alkalinity	Acidic
Aggregate form	Liquid concentrate
Composition	Citric acid
Biodegradable	•
Phosphate-free	•
Particularly effective at low temperatures	•
Highly effective, even at 40 °C	•
Short activation time	•
High-level detergency combined with sparing use	•
Density (at 20 °C) in g/cm ³	1.23
Based on 100 % sustainable raw materials	•
Degrades without any residues	•
Free from micro-plastics	•
Simple handling	•
Ideal additive for all detergents	•
Storage information	
Minimum storage temperature in °C	0
Maximum storage temperature in °C	40
Store in a dry location that is free of frost	•
Containers should always be re-sealed immediately once the required amount of product has been removed	•
Can be used after defrosting without loss of quality	•
Minimum shelf life in months	24
Dimensions and weight	
External dimensions, net height in mm	249
External dimensions, net width in mm	148
External dimensions, net depth in mm	193
External dimensions, gross height in mm	249
External dimensions, gross width in mm	148
External dimensions, gross depth in mm	193

ProCare Tex 30 C - 5 l

Neutralising agent, liquid concentrate, acidic, 5 l

for optimum protection of textiles thanks to reliable neutralisation.



EAN: 4002516776215 / Material number: 12495200 / Old Material Number: 59711101V1

Gross weight in kg	5.55 kg
--------------------	---------

Fill volume in ml	5000
-------------------	------

Compliance, test certificates, and symbols	
---	--

GHS/CLP compliant	•
-------------------	---