

## ProCare Tex 30 C - 5 l

Neutralising agent, liquid concentrate, acidic, 5 l  
for optimum protection of textiles thanks to reliable neutralisation.



- Neutralising agent based on citric acid
- For acidic pre-wash or for use after main wash
- Suitable for all wash processes
- Prevents yellow stains
- Stable, neutral pH values in the textiles

EAN: 4002516776215 / Material number: 12495200 / Old Material Number: 59711101V1

Product category	
Washing machines	•
Dispensing	
Suitable for automatic dispensing	•
Manual dispensing	•
Can be used with all degrees of water hardness	•
Dispensing recommendation	
pH value at 20 °C	2,5
How to use	
Launderettes	•
Hotels, catering	•
Cleaning sportswear	•
Cleaning outerwear	•
Cleaning workwear	•
Product properties	
Medium type	Neutralising agents
Alkalinity	Acidic
Aggregate form	Liquid concentrate
Composition	Citric acid
Biodegradable	•
Phosphate-free	•
Particularly effective at low temperatures	•
Highly effective, even at 40 °C	•
Short activation time	•
High-level detergency combined with sparing use	•
Density (at 20 °C) in g/cm <sup>3</sup>	1.23
Based on 100 % sustainable raw materials	•
Degrades without any residues	•
Free from micro-plastics	•
Simple handling	•
Ideal additive for all detergents	•
Storage information	
Minimum storage temperature in °C	0
Maximum storage temperature in °C	40
Store in a dry location that is free of frost	•
Containers should always be re-sealed immediately once the required amount of product has been removed	•
Can be used after defrosting without loss of quality	•
Minimum shelf life in months	24
Compliance, test certificates, and symbols	
GHS/CLP compliant	•