

A 315/1 Module

For increased capacity and reprocessing of lumened instruments.



- Increases the injector capacity from 22 to 44
- Can be used in all A 151 lower baskets available on the market
- Compatible with all standard dental & ENT accessories

EAN: 4002516881186 / Material number: 12531630 / Old Material Number: 69631501D

Product category • Thermal disinfectors • Washer-disinfectors, medical • Machine affiliation • PG 8581 • PG 8562 • PG 8582 CD • PG 8591 • PG 8592 • PWD 8682 • PWD 8682 CD • PWD 8682 CD • PWD 8692 •
Washer-disinfectors, medical • Machine affiliation • PG 8581 • PG 8562 • PG 8582 • PG 8582 CD • PG 8591 • PG 8592 • PWD 8682 • PWD 8682 CD • PWD 8682 CD •
Machine affiliation PG 8581 • PG 8562 • PG 8582 • PG 8582 CD • PG 8591 • PG 8592 • PWD 8682 • PWD 8682 CD • PWD 8682 CD •
PG 8581 PG 8562 PG 8582 PG 8582 PG 8582 CD PG 8591 PG 8592 PWD 8682 PWD 8682 C PWD 8682 CD PWD 8682 CD PWD 8682 CD PWD 8682 CD
PG 8562 PG 8582 PG 8582 CD PG 8591 PG 8592 PWD 8682 PWD 8682 CD
PG 8582 PG 8582 CD PG 8581 PG 8591 PG 8592 PWD 8682 PWD 8682 C PWD 8682 CD PWD 8682 CD PWD 8682 CD
PG 8582 CD PG 8591 PG 8592 PWD 8682 PWD 8682 C PWD 8682 CD PWD 8682 CD PWD 8682 CD
PG 8591 • • • • • • • • • • • • • • • • • • •
PG 8592
PWD 8682 • PWD 8682 C • PWD 8682 CD •
PWD 8682 C • PWD 8682 CD •
PWD 8682 CD •
DIV/D 8602
F WV D 0002
Product category
Insert for ENT instruments •
Insert/module for dental instruments •
Other inserts/other modules •
Product properties
Material Stainless steel, Plastic
Colour Stainless Steel, Black
Option to use a central filter •
Area of application
Reprocessing of ENT instruments
Reprocessing of dental instruments •
Reprocessing of lumened instruments •
Capacity
ENT instruments [number] 22
Transmission instruments [number] 22
Siphons [number] 22
Lumened instruments [number] 22
Dimensions and weight
External dimensions, gross height in mm 165
External dimensions, gross width in mm 505
External dimensions, gross depth in mm 550
Net weight in kg 2.90
Gross weight in kg 3.66
Companion pack
A 838 set of spanners [number]