

## TU HyClean Pure

HyClean Pure TU vacuum cleaner bag

Miele AirClean system delivers the best in filtration performance at 99.99 %.

EAN: 4002516829751 / Material number: 12557060 / Old Material Number: 41996682INT



- 80 % recycled material and ultimate performance – unrivalled on the market
- Exceptional service life thanks to patented 3D airflow guides
- HyClose automatic self-sealing collar for even greater hygiene
- Highest level of efficiency when vacuuming with a low wattage
- Contains: 4 vacuum cleaner bags, 1 AirClean exhaust filter
- Suitable for all Guard L1 and Guard S1 models