

APWD 338

Closure seals

with indicator to prevent unwanted opening after steam sterilisation.



- Indicator for steam sterilisation in acc. with EN ISO 11140-1 Class 1
- 100 units per pack

EAN: 4046963756512 / Material number: 12566280 /  
Old Material Number: 69633800D

Product category	
Small sterilisers	•
Machine affiliation	
PST 1710	•
PST 1720	•
PST 2210	•
PST 2220	•
Product category	
Product properties	
Material	Plastic
Colour	Orange
Area of application	
Reprocessing of ophthalmology instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of VET instruments	•
Reprocessing of podiatry instruments	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•