

MB1 250
 Measuring beaker, 250 ml
 for easy dispensing of ProCare products.



- Stackable measuring beaker with strong handle
- Printed with ml scale in increments of 5 ml
- For precise dispensing of liquid and powder products

EAN: 4002516805953 / Material number: 12594330 /
 Old Material Number: 59712001V1

Product category	
ProCare	•
Applications	
Easy dispensing of ProCare products	•
Dimensions and weight	
External dimensions, gross depth in mm	132.00
Dimensions and weight	
External dimensions, net height in mm	115
External dimensions, net width in mm	117
External dimensions, net depth in mm	77
External dimensions, gross height in mm	105
External dimensions, gross width in mm	94
External dimensions, gross depth in mm	132
Gross weight in kg	0