

## XKM 3200 WL LMD86

Communication module for PLW/PWD 868x or 869x generation machines for ExpertLine for connection to Miele MOVE, Segosoft or other software.

EAN: 4002516866527 / Material number: 12611070 / Old Material Number: 68822104D



- State-of-the-art connection to Miele MOVE
- Secure transmission via Ethernet
- REST interface with data provision in JSON format
- Flexible operation on local networks
- The last 30 process protocols are stored in the machine

Product category	
Washer-disinfectors, laboratory	•
Thermal disinfectors	•
Washer-disinfectors, medical	•
Washer-disinfectors, fire brigades and emergency services	•
Machine affiliation	
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
PLW 8683 [Safety]	•
PLW 8693 [Safety]	•
PWA 8672	•
PWD 8682	•
PWD 8682 C	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Process documentation/operating documentation/communication	•
Area of application	
Reprocessing of pipettes/phials	•
Reprocessing of laboratory glassware	•
Reprocessing of anaesthetic instruments	•
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Reprocessing of robotic instruments	•
Reprocessing of GYN instruments	•
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Reprocessing of theatre shoes	•
Reprocessing of basins	•
Reprocessing of kidney dishes	•
Reprocessing of baby bottles	•
Reprocessing of surgical instruments	•
Reprocessing of lumened instruments	•
Reprocessing of intramedullary nails/bone marrow drills	•
Reprocessing of trays	•
Reprocessing of mesh trays	•
Product properties	
Material	Plastic
Colour	Grey
Dimensions and weight	
External dimensions, net height in mm	40
External dimensions, net width in mm	90
External dimensions, net depth in mm	110
External dimensions, gross height in mm	105
External dimensions, gross width in mm	165
External dimensions, gross depth in mm	210

## XKM 3200 WL LMD86

Communication module for PLW/PWD 868x or 869x generation machines  
for ExpertLine for connection to Miele MOVE, Segosoft or other software.

EAN: 4002516866527 / Material number: 12611070 / Old Material Number: 68822104D



Net weight in kg	0.18
------------------	------

Gross weight in kg	0.46
--------------------	------