


APFD 201

Use

for the optimum loading of breathing masks.



- Suitable for holding 6 breathing masks securely
- Accessories: 1 small-part basket with 6 coloured compartments
- With separate lid for particularly delicate parts (6 compartments)

EAN: 4002516855248 / Material number: 12615430 / Old Material Number: 67020102D

Product category	
Fresh-water dishwashers	•
Machine affiliation	
PG 8063	•
PLW 8683 [Safety]	•
PLW 8693 [Safety]	•
Product category	
Insert for breathing masks	•
Area of application	
Reprocessing of breathing masks	•
Capacity	
Breathing masks [number]	6