

Guard L1

Cylinder vacuum cleaner

powerful | premium design | convenient cleaning with 2-in-1 nozzle



- ComfortFit: developed for the highest hygiene standards
- Durable: tested for the equivalent of 20 years** use
- Ideal power for every cleaning task: 4 selectable modes
- Individually selectable height and ergonomic comfort handle
- Premium design concept with stylish applications
- Range extra, fewer socket changes: 12 m operating radius
- 2-in-1 upholstery and crevice nozzle for easy cleaning without switching over
- Hygiene AirClean filter: effective filtration for clean room air

* Compared to SBB300-3 Parquet Twister parquet brush

EAN: 4002516833055 / Material number: 12652190 /
Old Material Number: 41UDF020GB

Construction type and design	
Colour	Nordic blue
Construction type	
Product type	SUDFO
User convenience	
Electronic suction power control	+/- footswitches
Illuminated power indicator	•
Park system for pauses during vacuuming	•
Park system	on both sides
Comfort cable rewind	•
Comfort handle 2.0	•
Mobility	
Castors with steel shaft, gentle on floors	•
3-point manoeuvrability	•
Operating radius in m	12
Cable length in m	8.5
Filter system	
Dustbag volume in l	3.5
Dustbag change indicator	•
Exhaust filter	AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic dustbag positioning	•
All-round furniture protection strip	•
Technical data	
Power	890
Volume in dB	75.0
Width in mm	500
Height in mm	239
Weight in kg	6.5
Accessories included	
Three-piece accessory set	Integrated
Upholstery nozzle	•
Crevice nozzle	•
Dusting brush with synthetic bristles	•
Floorhead	•
Vacuum cleaner bag voucher	•