

Guard L1 Flex
Cylinder vacuum cleaner
powerful | premium design | including flexible crevice nozzle



- 2-in-1 upholstery and crevice nozzle for easy cleaning without switching over
- Durable: tested for the equivalent of 20 years** use
- ComfortFit: developed for the highest hygiene standards
- Individually selectable height and ergonomic Comfort handle
- Ideal power for every cleaning task: 4 selectable modes
- Ideal for hard-to-reach places: flexible crevice nozzle
- AirClean filter: effective filtration for clean room air
- Premium design concept with stylish applications

* Compared to SBB300-3 Parquet Twister parquet brush

EAN: 4002516833062 / Material number: 12652200 /
Old Material Number: 41UDF030GB

Construction type and design	
Appliance colour	Sunset yellow
Construction type	
Cylinder vacuum cleaner with dustbag	•
Product brand	Guard L1
Motor type	PowerLine
Product type	SUDFO
User convenience	
Electronic suction power control	+/- footswitches
Illuminated power indicator	•
Suction tube	Suction tube X(HES)
Park system for pauses during vacuuming	•
Park system	on both sides
Comfort cable rewind	•
Comfort handle 2.0	•
Mobility	
Castors with steel shaft, gentle on floors	•
3-point manoeuvrability	•
Operating radius in m	12
Cable length in m	8.5
Filter system	
Vacuum cleaner bag type	HyClean Pure TU
Dustbag volume in l	3.5
Dustbag change indicator	•
Exhaust filter	AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic dustbag positioning	•
All-round furniture protection strip	•
Technical data	
Power	890
Volume in dB	74.0
Width in mm	500
Height in mm	239
Weight in kg	6.6
Accessories included	
Three-piece accessory set	Integrated
Upholstery nozzle	•
Crevice nozzle	•
Extra-long, flexible crevice nozzle SFD 20	•
Dusting brush	•
Floorhead SBD 365-3	•