

Guard L1 Allergy  
Cylinder vacuum cleaner  
powerful | HEPA AirClean filter | ideal for allergy sufferers



- 2-in-1 upholstery and crevice nozzle for easy cleaning without switching over
- Ideal power for every cleaning task: 4 selectable modes
- Premium design concept with stylish applications
- Individually selectable height and ergonomic comfort handle
- For effortless cleaning: 2-in-1 upholstery and crevice nozzle
- Ideal for allergy sufferers thanks to HEPA AirClean filter
- Durable: tested for the equivalent of 20 years\* use
- ComfortFit: developed for the highest hygiene standards

\* Compared to SBB300-3 Parquet Twister parquet brush

EAN: 4002516832706 / Material number: 12652230 /  
Old Material Number: 41UFF033GB

<b>Construction type</b>	
Product type	SUFFO
<b>Construction type and design</b>	
Colour	Brilliant white
<b>User convenience</b>	
Electronic suction power control	+/- footswitches
Illuminated power indicator	•
Park system for pauses during vacuuming	•
Park system	on both sides
Comfort cable rewind	•
Comfort handle 2.0	•
<b>Mobility</b>	
Castors with steel shaft, gentle on floors	•
3-point manoeuvrability	•
Operating radius in m	12
Cable length in m	8.5
<b>Filter system</b>	
Dustbag volume in l	3.5
Dustbag change indicator	•
Exhaust filter	HEPA AirClean filter
Motor protection filter	•
<b>Safety</b>	
Locking system for floorhead and handle	•
Automatic dustbag positioning	•
All-round furniture protection strip	•
<b>Technical data</b>	
Power	890
Volume in dB	74.0
Width in mm	500
Height in mm	239
Weight in kg	6.5
<b>Accessories included</b>	
Three-piece accessory set	Integrated
Upholstery nozzle	•
Crevice nozzle	•
Dusting brush with synthetic bristles	•
Floorhead	•
Vacuum cleaner bag voucher	•