

Guard L1 Cat & Dog Flex
Cylinder vacuum cleaner
powerful | attachment for removing pet hair | flexible crevice nozzle



- 2-in-1 upholstery and crevice nozzle for easy cleaning without switching over
- Active AirClean filter neutralises unpleasant odours
- Durable: tested for the equivalent of 20 years** use
- ComfortFit: developed for the highest hygiene standards
- For effortless cleaning: 2-in-1 upholstery and crevice nozzle
- Ideal power for every cleaning task: 4 selectable modes
- Ideal for hard-to-reach places: flexible crevice nozzle
- Individually selectable height and ergonomic comfort handle
- Premium design concept with stylish applications
- Turbobrush: removes stubborn dirt and animal hair

* Compared to SBB300-3 Parquet Twister parquet brush

EAN: 4002516833222 / Material number: 12652250 /
Old Material Number: 41UZFO35GB

Construction type	
Product type	SUZFO
Construction type and design	
Colour	Obsidian black
User convenience	
Electronic suction power control	+/- footswitches
Illuminated power indicator	•
Park system for pauses during vacuuming	•
Park system	on both sides
Comfort cable rewind	•
Comfort handle 2.0	•
Mobility	
Castors with steel shaft, gentle on floors	•
3-point manoeuvrability	•
Operating radius in m	12
Cable length in m	8.5
Filter system	
Dustbag volume in l	3.5
Dustbag change indicator	•
Exhaust filter	Active AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic dustbag positioning	•
All-round furniture protection strip	•
Technical data	
Power	890
Volume in dB	73.0
Width in mm	500
Height in mm	239
Weight in kg	6.5
Accessories included	
Three-piece accessory set	Integrated
Upholstery nozzle	•
Crevice nozzle	•
Extra-long, flexible crevice nozzle	•
Dusting brush with synthetic bristles	•
Floorhead	•
Vacuum cleaner bag voucher	•
Turbobrush STB 305-3	•