

Guard M1 Cat & Dog Flex Cylinder vacuum cleaner

powerful suction | attachment for removing pet hair | flexible crevice nozzle



- ComfortFit: developed for the highest hygiene standards
- Durable: tested for the equivalent of 20 years** use
- Ideal for hard-to-reach places: flexible crevice nozzle
- Active AirClean filter neutralises unpleasant odours
- Mode button: simply adjust the suction power to your needs
- Stylish design concept and illuminated LED icons
- Range extra, fewer socket changes: 11 m operating radius
- Turbobrush: removes stubborn dirt and animal hair

* Compared to SBB300-3 Parquet Twister parquet brush

EAN: 4002516833000 / Material number: 12652570 /
Old Material Number: 41VZF036GB

Construction type and design	
Colour	Obsidian black
Construction type	
Cylinder vacuum cleaner with dustbag	•
Product brand	Guard M1
Motor type	PowerLine
Product type	SVZF0
User convenience	
Electronic suction power control	MODE button
Ergonomic carrying handle	•
Ergonomic handle	•
Suction tube	HES mid-size telescopic tube
Park system for pauses during vacuuming	•
Park system	on both sides
Comfort cable rewind	•
Mobility	
Castors with steel shaft, gentle on floors	•
3-point manoeuvrability	•
Operating radius in m	11
Cable length in m	7.5
Filter system	
Vacuum cleaner bag type	HyClean Pure CO
Dustbag volume in l	4.5
Dustbag change indicator	•
Exhaust filter	Active AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic dustbag positioning	•
All-round furniture protection strip	•
Technical data	
Power	890
Volume in dB	77.0
Width in mm	496
Height in mm	227
Weight in kg	6.7
Accessories included	
Three-piece accessory set	Integrated
Upholstery nozzle	•
Crevice nozzle	•
Extra-long, flexible crevice nozzle	•
Dusting brush with synthetic bristles	•
Floorhead	•
Turbobrush STB 305-3	•