

Guard S1

Cylinder vacuum cleaner

high suction | light, compact | in modern trend colours | easy to adjust height



- ComfortFit: developed for the highest hygiene standards
- Ergonomic handle: for comfortable handling
- 99.99 %* dust retention capacity with AirClean filter
- Durable: tested for the equivalent of 20 years*** use
- 10 m operating radius: effortless cleaning of larger areas
- Easy handling thanks to lightweight and compact design
- Height-adjustable EasySlide telescopic tube

* In accordance with IEC 62885-2:2021

** Compared to SBB300-3 Parquet Twister parquet brush

EAN: 4002516833741 / Material number: 12652890 /

Old Material Number: 41TAF032GB

Construction type	
Product type	STAFO
Construction type and design	
Colour	Nordic blue
User convenience	
Electronic suction power control	4-stage rotary selector
Park system for pauses during vacuuming	•
Mobility	
Castors with steel shaft, gentle on floors	•
3-point manoeuvrability	•
Operating radius in m	10
Cable length in m	7.5
Filter system	
Dustbag volume in l	3.5
Dustbag change indicator	•
Exhaust filter	AirClean filter
Motor protection filter	•
Safety	
Locking system for floorhead and handle	•
Automatic dustbag positioning	•
All-round furniture protection strip	•
Technical data	
Power	890
Volume in dB	76.0
Width in mm	425
Height in mm	230
Weight in kg	5.7
Accessories included	
Three-piece accessory set	On VarioClip
Upholstery nozzle	•
Crevice nozzle	•
Dusting brush with synthetic bristles	•
Floorhead	•
Vacuum cleaner bag voucher	•