

U 403

Basket

for the optimum cleaning of cutlery.



- Close-meshed base
- Also suitable for sundae dishes
- Without cutlery inserts

EAN: 4002516854500 / Material number: 12654520 /
Old Material Number: 67640302D

Product category

Tank dishwashers

•

Machine affiliation

PG 8164

•

PTD 701

•

Product category

Basket for cutlery

•

Basket for saucepans/bowls, etc.

•

Product properties

Material

Wire, rilsan-coated

Colour

White

Dimensions and weight

External dimensions, gross height in mm

180

External dimensions, gross width in mm

520

External dimensions, gross depth in mm

610

Net weight in kg

1.85

Gross weight in kg

2.88