

U 404

Basket

For the optimum loading of cutlery.



- With 6 plastic cutlery holders
- Compartment size 106 x 106 mm

EAN: 4002516854647 / Material number: 12654560 /
Old Material Number: 67640402D

Product category

Tank dishwashers

•

Machine affiliation

PG 8164

•

PTD 701

•

Product category

Basket for cutlery

•

Product properties

Material

Wire, rilsan-coated, Plastic

Colour

Grey, White

Number of compartments

6

Compartment width in mm

106

Compartment length in mm

106

Dimensions and weight

External dimensions, gross height in mm

180

External dimensions, gross width in mm

520

External dimensions, gross depth in mm

610

Net weight in kg

1.62

Gross weight in kg

2.57