

U 407

Basket

For typical wash loads in bistros.



- For glasses, cups, plates, cake plates, saucers, etc.
- With different rail heights and distances between rails
- Quiver E 168 can also be used for cutlery

EAN: 4002516855255 / Material number: 12654670 / Old Material Number: 67640702D

Product category

Tank dishwashers

•

Machine affiliation

PG 8164

•

PTD 701

•

Product category

Basket for mixed tableware

•

Product properties

Material

Wire, rilsan-coated

Colour

White

Number of rows

5

Dimensions and weight

External dimensions, gross height in mm

180

External dimensions, gross width in mm

520

External dimensions, gross depth in mm

620

Net weight in kg

1.98

Gross weight in kg

2.93