


U 409

Basket

For typical wash loads in bistros.



- For glasses, cups, plates, cake plates, saucers, etc.
- Quiver E 165 can also be used for cutlery

EAN: 4002516854999 / Material number: 12654740 / Old Material Number: 67640902D

Product category	
Tank dishwashers	•
Machine affiliation	
PG 8164	•
PTD 701	•
Product category	
Basket for mixed tableware	•
Product properties	
Material	Plastic
Colour	Blue, Grey, White
Number of rows	3
Dimensions and weight	
External dimensions, gross height in mm	220
External dimensions, gross width in mm	550
External dimensions, gross depth in mm	550
Net weight in kg	1.47
Gross weight in kg	2.60