


E 167
Cutlery basket
For holding cutlery.



- Suitable for plate basket U 506
- 7 Compartments

EAN: 4002516854869 / Material number: 12663060 / Old Material Number: 67616702D

| Product category | |
|------------------------------|---|
| Fresh-water dishwashers | • |
| Tank dishwashers | • |
| Throughfeed tank dishwashers | • |
| Machine affiliation | |
| G 8066 | • |
| PG 8041 | • |
| PG 8043 | • |
| PG 8045 | • |
| PG 8055 | • |
| PG 8055 U | • |
| PG 8056 | • |
| PG 8056 U | • |
| PG 8057 TD | • |
| PG 8057 TD U | • |
| PG 8058 | • |
| PG 8058 U | • |
| PG 8059 | • |
| PG 8059 U | • |
| PG 8060 | • |
| PG 8061 | • |
| PG 8061 U | • |
| PG 8165 | • |
| PG 8166 | • |
| PG 8169 | • |
| PG 8172 | • |
| PG 8172 Eco | • |
| PTD 901 | • |
| PG 8096 | • |
| PG 8096 U | • |
| PG 8099 | • |
| PG 8099 U | • |
| PFD 400 | • |
| PFD 400 U | • |
| PFD 401 | • |
| PFD 401 U | • |
| PFD 402 | • |
| PFD 402 U | • |
| PFD 404 | • |
| PFD 404 U | • |
| PFD 405 | • |
| PFD 405 U | • |
| PFD 407 | • |
| PFD 407 U | • |
| PFD 408 U | • |
| PTD 702 | • |
| PTD 703 | • |
| PTD 704 | • |

E 167
Cutlery basket
For holding cutlery.



EAN: 4002516854869 / Material number: 12663060 / Old Material Number: 67616702D

| | |
|-----------------------------------------|---------|
| Product category | |
| Insert for cutlery | • |
| Product properties | |
| Material | Plastic |
| Colour | White |
| Number of compartments | 7 |
| Dimensions and weight | |
| External dimensions, gross height in mm | 95 |
| External dimensions, gross width in mm | 170 |
| External dimensions, gross depth in mm | 550 |
| Net weight in kg | 0.19 |
| Gross weight in kg | 0.34 |