

E 168
Cutlery basket
for the optimum loading of cutlery.



- Can be used in bistro basket U 407
- 4 compartments

EAN: 4002516854401 / Material number: 12663110 / Old Material Number: 67616802D

Product category	
Tank dishwashers	•
Machine affiliation	
PG 8164	•
PTD 701	•
Product category	
Insert for cutlery	•
Product properties	
Material	Wire, rilsan-coated
Colour	White
Number of compartments	4
Dimensions and weight	
External dimensions, gross height in mm	85
External dimensions, gross width in mm	150
External dimensions, gross depth in mm	275
Net weight in kg	0.49
Gross weight in kg	0.61