


A 104

Upper basket

For the optimum loading of instruments for MIS, arthroscopy and urology.



- Right side for holding 11 hollow instruments
- Additional Luer lock connectors for hollow instruments
- With spray arm: Mesh trays can be placed on the lower carrier

EAN: 4002516854807 / Material number: 12664810 / Old Material Number: 69610401D

Product category

Washer-disinfectors, medical

•

Machine affiliation

PG 8562

•

PG 8582

•

PG 8582 CD

•

PG 8592

•

Product category

Basket for MIS instruments

•

Injector basket

•

Upper basket

•

Product properties

Material

Stainless steel

Colour

Stainless Steel

Area of application

Reprocessing of MIS instruments

•

Capacity

DIN mesh trays [number]

1

Lumened instruments [number]

11

Dimensions and weight

External dimensions, net height in mm

332

External dimensions, net width in mm

528

External dimensions, net depth in mm

526

External dimensions, gross height in mm

420

External dimensions, gross width in mm

533

External dimensions, gross depth in mm

594

Net weight in kg

4.86

Gross weight in kg

6.89