


A 105/1

Upper basket

Height adjustable, without injector nozzles, for various instruments.



- Right-hand side: space for 22 injectors
- Left-hand side free for inserts
- Active internal drying of lumened instruments
- Active inner drying of dental handpieces
- A 800 tubular filter can be used in both injector manifolds

EAN: 4002516854814 / Material number: 12664820 / Old Material Number: 69610502D

<b>Product category</b>	
Thermal disinfectors	•
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Injector basket	•
Upper basket	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Option to use a central filter	•
<b>Area of application</b>	
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
<b>Capacity</b>	
Transmission instruments [number]	22
Lumened instruments [number]	22
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	220
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	540
Net weight in kg	4.87
Gross weight in kg	5.56
<b>Companion pack</b>	
A 838 set of spanners [number]	1