


A 105/1

Upper basket

Height adjustable, without injector nozzles, for various instruments.



- Right-hand side: space for 22 injectors
- Left-hand side free for inserts
- #Active internal drying of handpieces and contra-angle handpieces
- Clearances H 225 +/- 30, W 360, D 445 mm
- A 800 tubular filter can be used in both injector manifolds

EAN: 4002516854814 / Material number: 12664820 / Old Material Number: 69610502D

Product category	
Thermal disinfectors	•
Cleaning machine, small, dental	•
Washer-disinfectors, medical	•
Machine affiliation	
PG 8581	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8591	•
PG 8592	•
PWA 8672	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Injector basket	•
Upper basket	•
Area of application	
Reprocessing of ENT instruments	•
Reprocessing of dental instruments	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Option to use a central filter	•
Capacity	
Transmission instruments [number]	22
Lumened instruments [number]	22
Dimensions and weight	
External dimensions, net height in mm	202
External dimensions, net width in mm	528
External dimensions, net depth in mm	528
External dimensions, gross height in mm	220
External dimensions, gross width in mm	540
External dimensions, gross depth in mm	540
Net weight in kg	4.88
Gross weight in kg	5.57
Companion pack	
A 838 set of spanners [number]	1