

A 860  
Holder rack



EAN: 4002516855118 / Material number: 12665410 / Old Material Number: 69686001D

<b>Product category</b>	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
<b>Machine affiliation</b>	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 7111	•
<b>Product category</b>	
Other inserts/other modules	•
<b>Product properties</b>	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Number of compartments	8
Compartment width in mm	114
Compartment length in mm	114
Part of the EasyLoad system	•
<b>Area of application</b>	
Reprocessing of laboratory glassware	•
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	170
External dimensions, gross width in mm	255
External dimensions, gross depth in mm	490
Net weight in kg	0.48
Gross weight in kg	0.88