

A 150

Lower basket

For the optimum loading of injector modules or inserts.



- Lower basket with two module connections
- for holding up to 2 injector modules or inserts
- Automatic closing of the connections when not in use
- Clearances: width 494, length 440 mm

EAN: 4002516854449 / Material number: 12666330 / Old Material Number: 69615001D

Product category	
Washer-disinfectors, laboratory	•
Machine affiliation	
PLW 8636 [LAB / LAB MON]	•
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 8604	•
PLW 8683	•
PLW 8683 CD	•
PLW 8693	•
Product category	
Basket with module connection	•
Lower basket	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Area of application	
Reprocessing of pipettes/phials	•
Reprocessing of laboratory glassware	•
Reprocessing of butyrometers	•
Reprocessing of measuring cylinders	•
Dimensions and weight	
External dimensions, gross height in mm	175
External dimensions, gross width in mm	545
External dimensions, gross depth in mm	565
Net weight in kg	2.63
Gross weight in kg	3.29