


A 201

Mobile unit

For the optimum loading of anaesthetic sets.



- 8 nozzles with spring holders for breathing tubes up to 1.5 m
- Storage on spiral
- 15 injector nozzles for intubation material

EAN: 4002516854005 / Material number: 12667090 / Old Material Number: 69620101D

|  |                          |
|--|--------------------------|
| <b>Product category</b>  |                          |
| Washer-disinfectors, medical                                     | •                        |
| <b>Machine affiliation</b>                                       |                          |
| PG 8562  | •                        |
| PG 8582  | •                        |
| PG 8582 CD   | •                        |
| PG 8592  | •                        |
| PWD 8682   | •                        |
| PWD 8682 CD  | •                        |
| PWD 8692   | •                        |
| <b>Product category</b>  |                          |
| Mobile unit for anaesthetic instruments                          | •                        |
| <b>Product properties</b>  |                          |
| Material   | Stainless steel, Plastic |
| Colour   | Stainless Steel, Black   |
| <b>Capacity</b>  |                          |
| AN sets [number]   | 3                        |
| <b>Dimensions and weight</b>                                     |                          |
| External dimensions, gross height in mm                          | 790                      |
| External dimensions, gross width in mm                           | 715                      |
| External dimensions, gross depth in mm                           | 860                      |
| Net weight in kg   | 8.97                     |
| Gross weight in kg   | 10.27                    |
| <b>Companion pack</b>  |                          |
| E 430/1 1/3 mesh tray [number]                                   | 1                        |
| E 431 injector nozzle for concertina hose, Ø 8 x 193 mm [number] | 1                        |
| E 432 holder for concertina breathing tubes [number]             | 1                        |
| E 433 holder for silicone breathing tubes [number]               | 2                        |
| E 434 holder for paediatric breathing tubes [number]             | 1                        |
| E 466 injector nozzle 8 x 333 mm, for breathing bags [number]    | 6                        |
| E 496 injector nozzle 4 x 120 mm [number]                        | 15                       |