


A 9/1

Perforated tray inlay

For the gentle cleaning of delicate items by interrupting the spray jet.



- 7 x 7 mm perforation, 3 mm slat

EAN: 4002516854029 / Material number: 12667330 /
Old Material Number: 69500903D

Product category

Large-capacity washer-disinfectors, laboratory

•

Machine affiliation

PG 8527

•

Product category

Underlay

•

Area of application

Reprocessing of laboratory glassware

•

Product properties

Material

Stainless steel

Colour

Stainless Steel

Dimensions and weight

External dimensions, net height in mm

1

External dimensions, net width in mm

773

External dimensions, net depth in mm

573

External dimensions, gross height in mm

10

External dimensions, gross width in mm

580

External dimensions, gross depth in mm

740

Net weight in kg

1.69

Gross weight in kg

2.08