

E 374

Insert

for optimum loading of ENT instruments such as inhalation connections.



- 24 comp. approx. 45 x 45 mm, 27 compartments approx. 12 x 12 mm
- Mesh width on base: 1.7 mm
- For use in upper and lower basket

EAN: 4002516854050 / Material number: 12667400 / Old Material Number: 69537402D

<b>Product category</b>	
Washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Insert for ENT instruments	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	51
<b>Area of application</b>	
Reprocessing of ENT instruments	•
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	88
External dimensions, gross width in mm	225
External dimensions, gross depth in mm	468
Net weight in kg	1.14
Gross weight in kg	1.34