

A 203

Mobile unit

for the optimum loading of MIS instruments.



- Suitable for holding and for E 451, E 457, E 460, E 473, E 444
- Loading on 2 levels

EAN: 4002516853985 / Material number: 12667770 / Old Material Number: 69620301D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Mobile unit for MIS instruments	•
Area of application	
Reprocessing of MIS instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Dimensions and weight	
External dimensions, net height in mm	520
External dimensions, net width in mm	529
External dimensions, net depth in mm	542
External dimensions, gross height in mm	540
External dimensions, gross width in mm	560
External dimensions, gross depth in mm	560
Net weight in kg	9.16
Gross weight in kg	11.00
Companion pack	
E 336 MIBO irrigation sleeve, Ø 11 mm, length 121 mm [number]	3
E 362 blind stopper [number]	2
E 442 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 6 mm [number]	15
E 443 irrigation sleeve, Ø 11 mm, length 121 mm, with cap Ø 10 mm [number]	5
E 444 drum for fibre optic cables/suction hoses [number]	1
E 445 – 12 covers for E 442 irrigation sleeve, opening with Ø 6 mm [number]	1
E 446 – 12 covers for E 443 irrigation sleeve, opening with Ø 10 mm [number]	1
E 447 female adapter, for male Luer Lock [number]	3
E 448 silicone hose with female Luer Lock adapter [number]	6
E 449 male adapter, for female Luer Lock [number]	5
E 451 insert 1/6 mesh tray [number]	1
E 452 injector nozzle 2.5 x 60 mm [number]	3
E 453 injector nozzle 4 x 110 mm with retaining clip [number]	8
E 454 injector nozzle for trocar sleeves 10–15 mm [number]	6

A 203
Mobile unit
for the optimum loading of MIS instruments.

EAN: 4002516853985 / Material number: 12667770 / Old Material Number: 69620301D



E 456 spring stay for MIS instruments [number]

4

E 464 holder for injector nozzle E 454 [number]

3