

A 301/4

Injector module

for optimum loading of laboratory glassware – can be equipped with 18 nozzles.

EAN: 4002516854739 / Material number: 12668790 / Old Material Number: 69630106D



- Suitable for lab. glassware incl. Erlenmeyer and round flasks.
- Module w/o injector nozzles – equipped based on individual needs
- For up to 18 injector nozzles
- Part of the EasyLoad system
- Nozzle spacing 78 mm across the width, 78 mm across the depth
- Particularly suitable for A 842 and A 843 (nozzle diameter 4 mm)

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine affiliation	
PG 8504	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
PLW 7111	•
Product category	
Insert/module for narrow-necked glassware	•
Insert for various flasks	•
Product properties	
Material	Stainless steel, Plastic
Colour	Stainless Steel, Black
Part of the EasyLoad system	•
Area of application	
Reprocessing of laboratory glassware	•
Capacity	
Narrow-necked glassware [number]	18
Dimensions and weight	
External dimensions, gross height in mm	100
External dimensions, gross width in mm	205
External dimensions, gross depth in mm	487
Net weight in kg	1.68
Gross weight in kg	1.89