

E 416
Insert
for optimum loading of 6 single and two-piece specula.



- Suitable for 6 single and two-part specula
- 7 holders 40 mm apart
- For use in upper (only with single-part specula) and lower basket

EAN: 4002516854135 / Material number: 12669750 / Old Material Number: 69541602D

Product category	
Washer-disinfectors, medical	•
Machine affiliation	
G 7831	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8534	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Insert for gynaecological instruments	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Area of application	
Reprocessing of GYN instruments	•
Capacity	
GYN specula [number]	6
Dimensions and weight	
External dimensions, net height in mm	157
External dimensions, net width in mm	178
External dimensions, net depth in mm	279
External dimensions, gross height in mm	210
External dimensions, gross width in mm	200
External dimensions, gross depth in mm	365
Net weight in kg	0.69
Gross weight in kg	0.97