

E 458

Insert

for the optimum loading of screw-cap teats.



- Suitable for 36 screw connection teats
- 36 compartments at 22 x 22
- equipped with a flip-up and lockable lid

EAN: 4002516854265 / Material number: 12669800 /
Old Material Number: 69545802D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7826	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Other inserts	•
Area of application	
Reprocessing of baby bottles	•
Product properties	
Material	Stainless steel
Colour	Stainless Steel
Number of compartments	36
Dimensions and weight	
External dimensions, net height in mm	63
External dimensions, net width in mm	215
External dimensions, net depth in mm	445
External dimensions, gross height in mm	90
External dimensions, gross width in mm	230
External dimensions, gross depth in mm	465
Net weight in kg	1.50
Gross weight in kg	1.77