


A 605  
Module  
for optimum reprocessing of vials and test tubes, etc.



- Equipped with 98 injector nozzles
- Max. H 160 mm, max. Ø 20 mm
- For use in A 503 mobile unit

EAN: 4002516854074 / Material number: 12669870 / Old Material Number: 69660501D

<b>Product category</b>	
Large-chamber laboratory washers	•
<b>Machine affiliation</b>	
PLW 8615	•
PLW 8616	•
PLW 8617	•
<b>Product category</b>	
Insert/module for test tubes	•
<b>Area of application</b>	
Reprocessing of pipettes/phials	•
Reprocessing of test tubes	•
Reprocessing of laboratory glassware	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Capacity</b>	
Test tubes [number]	98
Phials [number]	98
<b>Dimensions and weight</b>	
External dimensions, net height in mm	158
External dimensions, net width in mm	254
External dimensions, net depth in mm	604
External dimensions, gross height in mm	190
External dimensions, gross width in mm	288
External dimensions, gross depth in mm	620
Net weight in kg	2.85
Gross weight in kg	3.31