

E 466

Injector nozzle

for the optimum loading of breathing bags 8 x 333 mm.



EAN: 4002516854920 / Material number: 12669890 /  
Old Material Number: 69746602D

<b>Product category</b>	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
<b>Machine affiliation</b>	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
<b>Product category</b>	
Other accessories	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless Steel
<b>Area of application</b>	
Reprocessing of anaesthetic instruments	•
<b>Dimensions and weight</b>	
External dimensions, gross height in mm	12
External dimensions, gross width in mm	150
External dimensions, gross depth in mm	200
Net weight in kg	0.85
Gross weight in kg	1.00