

E 479
Holder
for micro instruments up to 4 mm Ø.



- For use in trays with 5 mm mesh
- 50 items per bag

EAN: 4002516854067 / Material number: 12669910 / Old Material Number: 69747902D

Product category	
Washer-disinfectors, medical	•
Large-capacity washer-disinfectors, medical	•
Machine affiliation	
G 7823	•
G 7824	•
G 7825	•
G 7826	•
PG 8527	•
PG 8528	•
PG 8535	•
PG 8536	•
PG 8562	•
PG 8582	•
PG 8582 CD	•
PG 8592	•
PWD 8682	•
PWD 8682 CD	•
PWD 8692	•
Product category	
Holder	•
Product properties	
Material	Plastic
Colour	Transparent
Area of application	
Reprocessing of MIS instruments	•
Reprocessing of ophthalmology instruments	•
Dimensions and weight	
External dimensions, gross height in mm	12
External dimensions, gross width in mm	100
External dimensions, gross depth in mm	150
Net weight in kg	0.07
Gross weight in kg	0.08